






# Clik-FIX® All-in-One Epidural / PNB Catheter Securement Device

## REF ECS-2000B

### DESCRIPTION

The B. Braun Medical Clik-FIX All-in-One Epidural / PNB Catheter Securement Device is an integrated securement device (ISD) designed to protect the insertion site and secure the catheter. The Clik-FIX Epidural / PNB device is compatible with 18 - 24 Ga. epidural and peripheral nerve block catheters.

### INDICATIONS

The Clik-FIX All-in-One Epidural / PNB Catheter Securement Device is indicated for use with 18 - 24 Ga. epidural and peripheral nerve block catheters.

### CONTRAINDICATIONS

Do not use the device on compromised skin (i.e., open wounds, burns, lesions, rash, or other irritation and/or patients with known tape or adhesive allergies or hypersensitive skin).

### ADVERSE REACTIONS

Skin irritation.

### WARNING AND PRECAUTIONS

- Follow facility protocol when placing and removing the device. Observe universal precautions.
- Do not use a Clik-FIX device on patients who are allergic to adhesives and on non-adherent skin.

- Do not apply alcohol, acetone, or other liquid agent directly on the Clik-FIX device when in use, as it can weaken the device.
- During application and removal of the Clik-FIX device, take care to minimize manipulation and movement of the catheter.
- Prepare the skin surface according to facility protocol. Remove oil and moisture from application area and allow the area to completely dry.
- Routinely inspect the Clik-FIX device, insertion site, and catheter position according to facility protocol and best practice guidelines.
- Replace the Clik-FIX device at least every 3 days or per facility protocol.

### INSTRUCTIONS FOR USE

#### 1. PREP

Establish catheter access according to facility protocol or remove previous dressing from catheter. Prepare and cleanse the insertion site and surrounding securement area according to facility protocol. Allow the area to completely dry.

#### 2. ALIGN AND PLACE

While stabilizing the catheter, align the device so the cut-out opening is facing away from the catheter trajectory and the insertion site is in the center of the cutout. Remove the center split release paper from the underside of the device one side at a time (Figure 2). Lightly press to adhere to the skin.

#### 3. SECURE

Carefully guide the catheter into the center securement channel (Figure 3a), lightly pressing it into place. Loop the catheter to the side and guide it down into the other channel (Figure 3b). Continue to guide the catheter around the bottom of the device and up through the remaining channel (Figure 3c). NOTE: Do not pass the catheter through each channel more than once.

#### 4. PROTECT

Using the paper tab on the bordered transparent dressing, peel off the release paper and align the dressing to be centered over the securement device (Figure 4a). Ensure that the insertion site and securement device are completely covered. Press lightly to adhere to the securement pad and skin (Figure 4b). Date and time per facility protocol.

### REMOVAL

#### 1. REMOVE TRANSPARENT DRESSING

While following facility protocol and stabilizing the catheter and skin, carefully remove the transparent dressing.

#### 2. REMOVE CATHETER FROM SECUREMENT CHANNELS

While stabilizing the catheter, carefully remove the catheter from the channels of the securement base pad.

#### 3. REMOVE DEVICE

While stabilizing the catheter, slowly peel back the edge of the securement base pad. If resistance is met, use an alcohol prep pad at the interface between the pad and skin surface to weaken the adhesive.

**STERILE EO**

Sterilized using ethylene oxide

Do not use if package is damaged

Do not resterilize

Single use

Consult Electronic Instructions for Use  
www.bbraunusa.com

Rx only

A4809212-1 REV. 4/24

Patent: [www.bbraunusa.com/patents](http://www.bbraunusa.com/patents)

Not made with natural rubber latex or phthalates.

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B. Braun Medical Inc.